BioLas® 6

Softlaser

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GML German Medical Laser GmbH

Contents

Intended Use.....	23
Preface.....	24
Safety precautions	25
Product description	26
Suggestions for use.....	27
BioLas® 6 start-up	28
Replacing the batteries	28
Surface treatment.....	29
What you need to know before using laser acupuncture ...	30
Laser acupuncture	30
Cleaning the laser tip	33
Technical data	33
Shipping and storage operating instructions.....	34
Operating conditions	34
Checks before each use	34
Warnings.....	34
Trouble shooting	35
Explanation of symbols used	35
Standards and regulations	36
Contraindications.....	36
Warranty and service life	36
Disposal	37
Recommended reading.....	38
Address of the manufacturer and service center	38
Manufacturer's declaration.....	39

Intended Use

The BioLas® 6 Softlaser is used to treat the indications given in these instructions for use.

Please also refer to the instructions for use and safety information, for use in accordance with regulations.

Before Use

Check the scope of delivery

Package content:

- 1 piece BioLas® 6 Softlaser
- 4 pieces batteries Typ AAA LR03
- 1 piece Operating Instructions
- 1 piece Guidebook for Softlaser-Therapy
- 1 piece Hard Case

Special accessories:

Eartip for Ear-Acupuncture
(optionally available, not included)

Preface

Congratulations on purchasing this high quality Softlaser.

Softlaser is an instrument which was mainly constructed to be easy and safe for people to use.

The radiation that comes from the Softlaser is controlled so that it neither generates heat nor cause human tissue to change.

The wave length of the laser beam is determined in such a way that the skin and the connective tissue are effectively and protectively penetrated

LASERRADIATION
DO NOT STARE INTO BEAM
EVEN WITH OPTICAL
INSTRUMENTS
LASER CLASS 2M



Safety Precaution


Please read carefully the operating instructions and safety precautions instructions for use.

- Warning! Danger for your eyesight
- Laser radiation can cause retinal damage even in the lower power range. Never point any laser at the eyes and never look directly into any laser beam, including through optical instruments such as glasses, magnifying glass, etc.
- Keep BioLas® 6 Softlaser out of reach of children. Children cannot be aware of the risk posed by a laser device.
- Always pass on your BioLas® 6 Softlaser to third parties only with the operating instructions.
- Do not modify the medical device.
- Do not use the device if there is any sign of damage
- Wounds and bacterial infections may only be treated with surface laser beams at a distance of 1 to 2 cm.
- Do not use the device in oxygen-rich surroundings or in combination with flammable gases.

Product Description

Unscrew the lasertip

Skintype

dry skin | normal skin | moist skin

Note:

To clean the laser tip (laser outlet), undo the tip in an anti-clockwise direction and clean the laser outlet with a brush, cotton swab or a soft cloth. The laser outlet opening must always be free. Then screw the laser tip back on to the appliance.

Attention! To avoid damage, turn on the tip slightly by hand.

Suggestions for use

- Acne /spots
- Eczema
- Arthritis
- Headaches
- Abscess (suppurative focus)
- Sinusitis
- Fever blisters
- Herpes Simplex
- Scars
- Wounds
- Muscle tension
- Pustules
- Chapped skin
- Gingivitis (Inflamed gums)
- Pulled ligaments-/ muscles
- Gout
- 1 and 2 grade burns
- Epicondylitis (tennis elbow)
- Tendinitis

And much more

BioLas® 6 corresponds to laser protection class 2M and has a visible laser beam and a variable laser beam cone.

You do not need to wear safety glasses or other protection.

The eye's closing reflex provides it with sufficient protection. Damage can only occur if you intentionally look into the beam for a long time. This should always be avoided.

BioLas® 6 Start-up

Inserting or changing the batteries (See “Product Description”, →p.26)

Unscrew the battery-holder screw until the battery-holder (9) can be removed from the housing. Insert the batteries in the holder in accordance with the correct polarity.

Push the battery-holder (9) back in the housing and carefully and firmly screw on the battery-holder screw(10).

Attention, only slightly tighten by hand

Replacing the batteries /Used batteries

If the control lamp (2) turns red, you must not replace the batteries. However, you should keep spare batteries. If the control lamp (2) is blinking red and the BioLas® 6 turns off, the batteries are used and should be replaced.

Remove the used batteries and replace them with four new 1.5 V alkaline manganese batteries, type AAA.

Attention!

Use only branded batteries protected against leakage! “Cheap” batteries can leak and damage the electronics!

If BioLas® 6 is not used for a long time, remove the batteries from the battery-holder.

Change used batteries immediately and dispose of them in an environmentally friendly way (See notice for battery disposal→p.37)

Surface treatment

Operation

(See „Product description“, → p.26.)

Switch the appliance on using the mains switch (1).

The control lamp (2) will light up green. The

BioLas® 6 is now ready for operation.

Now hold the BioLas® 6 in your hand as you would hold a pencil (see diagram) and point the laser tip (6) at the area to be treated, at a distance of 1 to 2cm.

Do not touch the point that requires treatment with the laser tip if you are treating surface damage to the skin (e.g. herpes, wounds etc).

Without Timer

Press the laser operation button (8) once with your index finger (if you are right handed) or thumb (if you are left handed). The laser is turned on. After the prescribed radiation time, press the laser operation button again. This will switch the laser off.

Note: in this mode there is no indication in the scale.

With Timer

Press and hold the laser-operating button . The LED-light scale shows now ascending the laser operation in minutes. 1 light = 1 min., 2 lights = 2 min. etc. Once you have reached the desired therapy time, release the laser-operating button. The laser turns on and the therapy begins. After each minute a light goes out. Upon with the expiration of the last light, the laser switches off automatically. You can turn off the laser during the operation at any time by pressing the laser operation-button again. The pre-programmed therapy time will be deleted. If you are using the laser for longer,

the laser will switch off automatically after 20 min. Please ensure that you switch off the appliance at the mains switch (1) after treatment, because electronic appliances use electricity constantly, even in standby mode, which could result in the batteries being empty within a few days.

What you need to know before using laser acupuncture

Because everyone is different, laser acupuncture will have a different effect depending on the individual. It is also impossible to give precise information on the treatment

Give ray treatment to the acupuncture points one after another by setting the laser tip to the treatment area and then start the laser as described below.

To find an acupuncture point at another person, you must keep up with the person's skin contact.

You will find some examples of acupuncture point areas in the GML Guidebook for Softlaser-Therapy or in your acupuncture books.

Laser acupuncture

Operating

(See „Product description“, → p.26.)

1. Switching on

Switch the appliance on using the mains switch (1). The control lamp (2) will light up green. The BioLas® 6 now ready for use and the point-detector is already activated.

2. Find out the skin type

Hold the BioLas® 6 in your hand like a pencil (see diagram). Touch the contact plate (3) with your thumb (if you are right handed) or your index finger (if you are left handed). If you do not touch the contact plate, there will be no display on the point-detector scale (5).

Now set the skin type adjustment switch (4) at the central position (2 drops = normal skin). Determine the skin type by placing the laser tip (6) vertically on the skin near any acupuncture point or an acupuncture point of your choosing (e.g. on your hand) with a light, constant pressure and gliding it around. The green point-detector scale (5) will now light up to various lengths (1 to 19 light points).

If no fields or only a few fields light up, the skin type is different.

If this is the case, set the skin type switch (4) at 1 drop = dry skin. If all of the fields of the green point-detector scale (5) light up constantly, you have moist skin. Set the skin type determination switch (4) at 3 drops = moist skin.

3. Point search and laser acupuncture

Now you can precisely determine the specific acupuncture point required for treatment. Proceed with the laser tip (6) as described above, over the area of the skin where the acupuncture point is located. The contact plate (3) must always be in contact with the thumb or index finger. The part of the skin where the point-detector scale (5)

shows the highest amplitude is the acupuncture point.

Now hold the BioLas® 6 directly above this point and briefly press the laser operation button (8) until the first led light on the point-detector scale (5) is shining.

If you release the laser operating button (8), the laser turns on and the acupuncture begins. The laser beam will switch off automatically after 1 min. and the point detector is reactivated, so that you can immediately find the next point and treat it.

During laser treatment, keep the BioLas® 6 as perpendicular to the surface of the skin as possible to get the best effect. The focal point of the laser is exactly located at the laser tip outlet.

By pressing the laser operation button (8), the point detector scale (5) will go out. The laser beam is visible. Therefore, the skin around the acupuncture point will be light up red during the acupuncture.

At the end of the laser acupuncture, please switch off at the mains switch (1). The control lamp (2) will go out.

Please ensure that you switch off the appliance at the mains switch (1) after treatment, because electronic appliances use electricity constantly, even in standby mode, which could result in the batteries being empty within a few days.

Changing the laser tip


To remove the laser tip undo the tip in an anti-clockwise direction. The screw on the ear tip (optional accessories) on to the appliance.

Warning! In order to avoid damage screw on the laser tip only by hand.

Cleaning the laser tip. (See product description“, → p.26.)

If the laser tip is pressed on the skin, it will become covered with skin particles or oil. This prevents the laser beam from working at full strength. Therefore clean the laser tip at regular intervals.

To clean the laser outlet, undo the laser tip in an anti-clockwise direction. Clean the laser tip with a brush, soft cloth or a cotton swab. Then screw the laser tip back onto the appliance.

Warning! In order to avoid damage only put on the laser tip by hand.

Clean the housing if necessary with a soft cloth. Use only mild cleaning agents.

Technical Data

Operating voltage: 6V direct current via 4 micro batteries type Micro LR03/AAA (1.5V)

Laser diode: visible red-light laser diode

Wavelength: Approx. 650nm

Beam divergence angle: Approx, 17 degrees.

Laser output: 5 mW, divergent continuous beam

Laser protection class: 2M

Medical appliance class: II a

Shipping and Storage Conditions

Ambient temperatures: from 0°C to +50°C above the melting point

Relative humidity: from 0% to 90% at +50°C, no condensation

Atmospheric pressure (barometric): from 500 kPa to 1060 kPa

Operating conditions

Temperature (recommended): from +18°C to 30°C above the melting point

Humidity (recommended): from 20% to 70% at +25°C, no condensation

Atmospheric pressure (barometric): from 800 kPa to 1060 kPa.

This corresponds to a working height of 0 to 2,000 meters above sea level

Checks before each use

Check the housing for damage

Laser tip free of dirt

Warnings

- Do not use if the housing is damaged
- Do not use if the battery holder is damaged
- Do not use in the shower or bath
- Do not open the housing

Troubleshooting

(See „Product description“, → p.26.)

The control lamp (2) is not shining → main switch is not in position I, the batteries are not inserted, discharged or not inserted correctly

Weak or no laserlight at the laser outlet → the outlet is dirty

No operation in spite of new batteries → the batteries are inserted incorrectly or defective

Explanation of Symbols Used

Symbol	Meaning	Location
	Laser warning label (yellow background) DO NOT LOOK DIRECTLY INTO THE LASER BEAM EVEN USING OPTICAL INSTRUMENTS LASER CLASS 1M IEC 60825-1	Laterally on the tip of the device
I O	ON/OFF switch button	On the top of the device
	Skintype settings for the Point-detector 1 drop = dry skin 2 drops = normal skin 3 drops = moist skin	Laterally on the device
	BF type operating unit	Laterally on the device
	Do not dispose of as household waste	Laterally on the device
	Timer	Laterally on the device

Standards and Regulations

BioLas ® 6 meets the requirements of EU Directive 93/42/EWG and the quality assurance requirements of EU Directive 93/42/EWG Annex VII.

The required harmonized standards are fulfilled

Contraindications

- Patients with pacemakers: no radiation treatment around the ribcage or pacemaker
- Pregnancy: no radiation treatment in the abdominal area / lower thorax
- Patients with tumours: no radiation treatment on the tumour or in case of suspicion of relapse or metastases
- No radiation treatment of Epiphysis (growth area of bones) and hollow bones of children
- No radiation treatment of the thyroids, ovaries and testicles
- No radiation treatment for epilepsy (oligoepilepsy)

Warranty and Service Life

BioLas ® 6 has a guarantee period of two years from the date of purchase (please present your purchase receipt).

During the warranty period, we will remedy any defects in the product free of charge, if the product is proven to be defective in workmanship or materials. The warranty does not cover high-wear parts. Moreover, the warranty is cancelled in case of improper use or non-observance of the operating instructions, use of aggressive cleaning agents, unauthorized intervention, or modifications of the device. When evaluating the medical product, its estimated service life is 10 years.

Disposal

Always observe your local regulations regarding disposal of waste available at your municipality or your dealer where you have purchased the product. All electrical appliances must be disposed of separately from regular household waste in designated areas.

Notice for battery disposal

In connection with the sale of batteries or with the delivery of devices containing batteries, we are obliged to draw your attention at the following: as consumer, you are legally obliged to return used batteries. You can return batteries after use to the seller or to the foreseen collection points (e.g. in local collection or commercially) free of charge. You can return the batteries by mail to the seller. The seller refunded to you in any case the postage for the return of the used batteries.

The symbols shown on the batteries have the following meaning:

The symbol of the crossed trashcan means that the battery should not be disposed of with household waste.

Pb = battery contains more than 0,004 percent plumb

Cd = battery contains more than 0,002 percent cadmium

Hg = battery contains more than 0,0005 percent mercury

Recommended Reading

You can find a wide range of specialized literature on laser therapy and laser acupuncture in German and English in bookshops.

Information about other products by GML for laser therapy is available in the Internet at:

www.gm-laser.de and
www.tinnimed.de

Address of the manufacturer and service center

GML German Medical Laser GmbH
Reutlinger Straße 14, D 75179 Pforzheim
Tel: + 49 (0) 7231 – 139569-0
Fax: + 49 (0) 7231 – 139569-50
Email: info@gm-laser.de

Manufacturer's Declaration

Electromagnetic compatibility (EMC)

Table 201

Guidelines and Manufacturer's Declaration – Electromagnetic emissions			
BioLas® 6 is designed to be operated in an environment as described below. The customer or user of BioLas® 6 should ensure that it is operated in such an environment.			
Emitted interference measurements	Compliance	Electromagnetic environment – Guidelines	
HF emissions pursuant to CISPR11	Group A	BioLas® 6 uses HF energy exclusively for its internal functions. Consequently, its HF emissions are very low, and it is unlikely that nearby electronic devices will be disturbed.	
HF emissions pursuant to CISPR11	Class B	BioLas® 6 is designed is suitable for use in facilities other than the residential area and facilities such as those directly connected to a public grid which also supplies buildings used for residential purposes.	

Table 202

Guidelines and Manufacturer's Declaration – Electromagnetic immunity			
BioLas® 6 is designed to be operated in an environment as described below. The customer or user of BioLas® 6 should ensure that it is operated in such an environment.			
Electromagnetic immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment / Guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Flooring should be made of wood or concrete, or be covered with ceramic tiles. If the floor is covered with synthetic material, the relative air humidity must be at least 30%.
Magnetic field at the supply frequency (50/60Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should conform to the typical values as they can be found in business and hospital environments.

Note: U_T = Public alternating current supply before the application of the test level.

Table 204

Guidelines and Manufacturer's Declaration – Electromagnetic immunity			
BioLas® 6 is designed to be operated in an environment as described below. The customer or user of BioLas® 6 should ensure that it is operated in such an environment.			
Electromagnetic immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment / Guidelines
			Portable and mobile radio equipment (incl. electric wiring) shall not be used any closer to BioLas® 6 than the recommended safe distance.
Conducted HF interference pursuant to IEC 61000-4-6	3 V 150 kHz - 80 MHz	3 V	$d = [3.5/V1] \text{ VP}$
Radiated HF interference pursuant to IEC 61000-4-6	3 V/m 80 MHz - 2.5 GHz	3 V/m	$d = [3.5/E1] \text{ VP}$ 80 MHz - 800 MHz $d = [7/E1] \text{ VP}$ 800 MHz - 2.5 GHz

			<p>P is taken as the nominal emission power of the transmitter in Watts (W) according to the specifications of the transmitter manufacturer and d as the recommended safety distance in meters (m). The field strength of a stationary radio transmitter is less^a at all frequencies than the compliance level according to a local investigation.^b Interference is possible in the vicinity of appliances bearing the left symbol:</p>
<p>Note 1. At 80 MHz and 800 MHz the higher value shall apply.</p>			
<p>Note 2. These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection by buildings, objects and people.</p>			
<p>^a The field strength of stationary transmitters, such as for example the base stations of radio telephones and public mobile telephone services, amateur radio stations, AM- and FM broadcasting- and television stations, is incapable, in theory of being precisely determined in advance. In order to be able to evaluate an electromagnetic environment after the installation of high frequency stationary transmitters, an examination of the location is to be recommended. When the measured field strength at the location of a BioLas[®] 6 exceeds the above compatibility level, then the normal operating conditions of BioLas[®] 6 at that particular location will have to be more closely examined. If unusual power emission characteristics are observed, then it may be necessary to take additional measures, e.g. the reorientation or relocation of BioLas[®] 6</p>			
<p>^b In the frequency range from 150 kHz to 80 MHz the field strength is less than [V1] V/m.</p>			

Table 206

Recommended safe distances between portable and mobile HF communication devices and BioLas[®] 6

BioLas[®] 6 is intended to be operated in an electromagnetic environment where the interference factors of high frequency emissions are under control. The customer or the user of BioLas[®] 6 can assist in avoiding electromagnetic interference by maintaining minimum distances between portable and mobile high frequency communication devices (transmitters) and BioLas[®] 6 as recommended below in accordance with the maximum power emissions of the communication device.

Nominal output of the transmitter (W)	Safe distance according to transmission frequency (m)		
	150 kHz - 80 MHz	80 MHz - 800 MHz	800 MHz – 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

In regard to transmitters whose nominal power output is not given in the above table, the distance can be calculated by using the equation where P represents the nominal power output of the transmitter in Watts W as per the specifications of the transmitter manufacturer.

NOTE 1. At 80 MHz and 800 MHz the higher frequency range shall apply.

NOTE 2. ISM frequency ranges (for industrial, scientific and medical purposes) between 150 kHz and 80 MHz vary from 6.765 MHz to 6.795 MHz; from 13.553 MHz to 13.567 MHz; from 26.957 MHz to 27.283 MHz and from 40.66 MHz to 40.70 MHz.

NOTE 3. Compliance levels in the ISM frequency ranges between 150 kHz and 2.5 GHz are determined to reduce the possibility of interference inducing by portable/ mobile communication equipment, if it is inadvertently introduced into the patient area. Therefore, to calculate the recommended safe distance in these frequency ranges an additional factor of 10/3 is used.

NOTE 4. These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection by buildings, objects and people.

GML German Medical Laser GmbH